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THROAT SELF-TEST TO DETECT THE OMICRON VARIANT.

Microbiologists recommend that you perform the self-test in the throat to detect the new Omicron variant of the coronavirus.

Our antigen self-tests have been specially tested to detect the Omicron variant of the coronavirus. It seems that Omicron is not only detected in the nose but also quickly and easily in the throat. If you only test in the nose, you can sometimes get a negative test result on your self-test while you are still carrying the coronavirus in the throat.

To get a better result, it is best to test in the throat as well. How do you do that?

To perform the test with some throat mucus, it is best not to eat or drink anything for 15-30 minutes. Then clear your throat and try to cough up some throat mucus and get it on your tongue. Then take the cotton swab and run a sweeping and twisting motion over the back of the tongue. After the throat swab, pass the same swab through each nostril one more time, following the directions in the package insert. You can then complete the test by running the cotton swab through the buffer and dripping the buffer with the nasal/pharyngeal mucus sample into the test according to the directions in the instructions for use.

January 2022.