2019-nCoV Ag Saliva Rapid Test Card
(Immunochromatographie)
Suitable for self-testing use

Catalogue Number:
0699C8X001 (1 Test/Kit) 0699C8X005 (5 Tests/Kit) 0699C8X020 (20 Tests/Kit)

INTENDED USE

The Test Card is a lateral flow immunoassay intended for the in vitro qualitative detection of 2019-nCoV antigen in human saliva specimens. The Test Card can be used for individuals with or without symptoms or other epidemiological reasons to suspect COVID-19 infection. A positive result indicates 2019-nCoV infection. Please quarantine yourself and contact a doctor. Additional testing is necessary. Positive results do not rule out bacterial infection or co-infection with other viruses. A negative result does not rule out COVID-19 infection. It does rule out 2019-nCoV infection. Negative results should be considered in the context of a patient's recent exposures, history, and the presence of clinical signs and symptoms consistent with COVID-19. Contact a doctor and confirm with a PCR test, if necessary.

For in vitro use only. Suitable for self-testing use.

SUMMARY AND EXPLANATION

The novel coronaviruses belong to the β genus, 2019-nCoV is an acute respiratory infectious disease. People are generally susceptible. Currently, the novel coronavirus are seen as a major threat to public health. Seroconverting infected people can also be an infectious source. Based on the current epidemiological investigation, the incubation period is 1 to 14 days, mostly 3 to 7 days. The main manifestations include fever, fatigue and dry cough. Nasal congestion, runny nose, sore throat, myalgia and diarrhea are found in a few cases.

PRINCIPLE OF THE TEST

This Test Card uses a lateral-flow sandwich to legally detect the antigen of novel coronavirus (2019-nCoV) in saliva samples. During detection, the gold labeled anti-2019-nCoV monoclonal antibody in the labeling band binds to the 2019-nCoV antigen in the sample to form a complex, and the reaction complex moves forward along the nitrocellulose membrane under the action of chromatography. It is captured by the anti-2019-nCoV monoclonal antibody precoated by the Test area (T). If the sample does not contain 2019-nCoV antigen, a red color line cannot be formed in the Test area (T). Regardless of whether the sample to be tested contains 2019-nCoV antigen, a red reaction line will always form in the quality Control area (C), if the test has been performed properly.

MATERIALS AND COMPONENTS

Materials provided with the test kits

REF Component	0699C8X001	0699C8X005	0699C8X020
Pouch/test card and desiccant	1	5	20

TEST PROCEDURE

Before test, please read the instructions carefully.

- Take the Test Card to equilibrate to room temperature.
- Unpack the aluminum foil bag. Place the Test Card horizontally on the table.
- Insert the saliva swab with collected sample into the Test Card holder and push down saliva swab. The bump at the end of the saliva swab must be into the hole of the Test Card holder.

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2. Limitations:

1. Das Test result of the test card should not be taken as a confirmed diagnosis, for clinical reference only. Judgement should be made along with RT-PCR results, clinical symptoms, epidemiological information, and laboratory data.

2. The kit are to be used for the qualitative detection of 2019-nCoV N-protein antigens from saliva specimen.

3. Test performance depends on the amount of virus antigen in the sample and may or may not correlate with viral culture results performed on the same sample.

4. The test card must stay at room temperature (15 - 30 °C) for 30 minutes before use, otherwise the results may be incorrect.

5. A false negative test result may occur if the level of antigen in a sample is below the detection limit of the test.

6. Failure to follow the Test Procedure may adversely affect test performance and/or invalidate the test result.

7. Reading time less than 10 minutes may lead to a false negative result; Reading time more than 15 minutes may lead to a false positive result.

8. Use appropriate precautions in the collection, handling, storage, and disposal of samples and used test contents.

9. Negative test results are not intended to rule in other viral or bacterial infections.

10. Negative results should be treated as presumptive and confirmed with a molecular assay.

11. Users should test specimens as quickly as possible after specimen collection.

12. If the sample volume is not enough, the test cannot be carried out successfully.

3. Interference Substances:
The test results do not be interfered with the substance at the following concentration:

No.	Specimen type	Conc.
1	HCoV-HKU1	10 ³ TCID ₅₀ /mL
2	Staphylococcus aureus	10 ³ CFU / mL
3	Streptococcus pyogenes	10 ³ CFU / mL
4	Measles virus	10 ³ TCID ₅₀ /mL
5	Paramyoxovirus parotitis	10 ³ TCID ₅₀ /mL
6	Adenovirus 3	10 ³ TCID ₅₀ /mL
7	Mycoplasma pneumoniae	10 ³ CFU / mL
8	Parainfluenza virus 2	10 ³ TCID ₅₀ /mL
9	Human Metapneumovirus (MPV)	10 ³ TCID ₅₀ /mL
10	Human coronavirus OC43	10 ³ TCID ₅₀ /mL
11	Human coronavirus 229E	10 ³ TCID ₅₀ /mL
12	Human coronavirus NL63	10 ³ TCID ₅₀ /mL
13	MERS-Coronavirus EMC/2012	10 ³ TCID ₅₀ /mL
14	Bordetella parapertussis	10 ³ CFU / mL
15	Influenza B (Victoria strain)	10 ³ TCID ₅₀ /mL
16	Influenza B (Y strain)	10 ³ TCID ₅₀ /mL
17	Influenza A (H1N1 2009)	10 ³ TCID ₅₀ /mL
18	Influenza A (H3N2)	10 ³ TCID ₅₀ /mL
19	Avian influenza virus (H7N9)	10 ³ TCID ₅₀ /mL
20	Avian influenza virus (H5N1)	10 ³ TCID ₅₀ /mL
21	Epstein-Barr virus	10 ³ TCID ₅₀ /mL
22	Enterovirus CA16	10 ³ TCID ₅₀ /mL
23	Rhinovirus	10 ³ TCID ₅₀ /mL
24	Respiratory syncytial virus	10 ³ TCID ₅₀ /mL
25	Streptococcus pneumoniae	10 ³ CFU / mL
26	Candida albicans	10 ³ CFU / mL
27	Chlamydia pneumoniae	10 ³ CFU / mL
28	Bordetella pertussis	10 ³ CFU / mL
29	Pneumocystis jirovecii	10 ³ CFU / mL
30	Mycobacterium tuberculosis	10 ³ CFU / mL
31	Legionella pneumophila	10 ³ CFU / mL

4. Interference Substances:
The test results do not be interfered with the substance at the following concentration:

No.	Interference substances	Conc.
1	White Blood	4%
2	Ibuprofen	1mg / mL
3	Tetracycline	3µg / mL
4	Chloramphenicol	3µg / mL
5	Erythromycin	3µg / mL
6	Tobramycin	5%
7	Throat spray (Menthol)	15%
8	Mupirocin	10mg/mL
9	Throat lozenges (Menthol)	1.5mg/mL
10	Tamiflu (Oseletamivir)	5mg/mL
11	Naphthoxoline hydrochloride nasal drops	15%
12	Mucin	0.50%
13	Fisherman's Friend	1.5mg/mL
14	Compound Benzocaine Gel	1.5mg/mL
15	Cromoglycate	15%
16	Sinex (Phenylephrine Hydrochloride)	15%
17	Afrin (Oxymetazoline)	15%
18	Fluticasone propionate spray	15%

5. Precision:

a. Test 100% indicates of negative and positive reference controls. The negative agreement and the positive agreement were 100%.

b. Test three different lots including positive and negative reference controls. The negative agreement and the positive agreement were 100%.

6. Hook Effect:
The Test Card was tested up to 1.6 × 10³ TCID₅₀/ml of heat-inactivated 2019-nCoV strain and no hook-effect was observed.

7. Limitations:

1. This product can only perform qualitative detection of 2019-nCoV antigen.

2. Positive Result:
If both C and T lines are visible, the test result is positive and valid.

3. Negative Result:
If a colored line is visible in the control area, and no colored line appears in the test area, the result is negative and valid.

4. Invalid Result:
The test result is invalid if a colored line does not form in the control area. The sample must be re-tested, using a new test card.

8. Performance Characteristic:

1. Clinical Verification:
The performance of Test Card was established with 793 saliva sample collected from symptomatic and asymptomatic patients.

2. Comparative RT-PCR Test Result:

	Comparative RT-PCR Test Result		
	Positive (+)	Negative (-)	Total
Detected Positive	258	6	264
Detected Negative	22	507	529
Total	280	513	793

3. Sensitivity:
92.14%, 95% CI (88.39, 94.75)

4. Specificity:
98.83%, 95% CI (97.47, 99.46)

5. Accuracy:
96.47%, 95% CI (94.94, 97.55)

9. Limit of Detection:
The experimental results show that for the virus culture concentration above 100 TCID₅₀/mL, the positive rate of detection is greater than or equal to 95%. So, the limit of detection of the Test Card is 100 TCID₅₀/mL.

10. Cross-reactivity:
Cross-reactivity of the test card was evaluated. The results showed no cross reactivity with the following specimen:

Action	Test result is positive:	Test result is negative:	Test result is invalid:
If the test result is positive:	- Note that you are currently suspected of COVID-19 infection. - Contact a doctor or the local health office immediately. Any medical decision should not be made without consulting a doctor or the local health office. - Comply with local self-isolation guidelines. - Do a PCR test for confirmation.	- Continue to comply with all local applicable rules and protective measures. - Be aware that even if the test result is negative, an infection may occur. - In case of doubt, repeat the test after 1-2 days, as the 2019-nCoV may not be detected accurately in all phases of an infection. - Do a PCR test for confirmation, if necessary.	- Possibility caused by incorrect operation. - Repeat a test, using a new test cassette. - If test result is still invalid, contact the distributor or the store where you bought the product, with the lot number.

11. Inadequate sample collection, processing, storage and transport may yield a false positive result or a false negative result.

12. To obtain accurate results, do not use visually cloudy or overly viscous samples.

13. Do not use the test card if it is damaged or expired.

14. Keep out of the reach of children.

15. The device should be operated by adult person.

16. Wear a safety mask or other face covering when collecting specimen from child or another individual.

17. Use of Nitrile, Latex (or equivalent) gloves is recommended when handling samples.

18. Wash hands thoroughly after handling.

19. Discard the diagnostic kit with the infectious waste.

20. Use local laws for disposal of the diagnostic kit.

21. Do not use the test card if it is damaged.

22. Do not use the test card if it is beyond its use-by date.

23. Do not use the test card if it is stored at temperatures higher than 30 °C.

24. Do not use the test card if it is stored at temperatures lower than 2 °C.

25. Do not use the test card if it is stored in direct sunlight.

26. Do not use the test card if it is stored near heat sources.

27. Do not use the test card if it is stored near strong magnetic fields.

28. Do not use the test card if it is stored near strong electromagnetic fields.

29. Do not use the test card if it is stored near strong electric fields.

30. Do not use the test card if it is stored near strong magnetic fields.

31. Do not use the test card if it is stored near strong electromagnetic fields.

32. Do not use the test card if it is stored near strong electric fields.

33. Do not use the test card if it is stored near strong magnetic fields.

34. Do not use the test card if it is stored near strong electromagnetic fields.

35. Do not use the test card if it is stored near strong electric fields.

36. Do not use the test card if it is stored near strong magnetic fields.

37. Do not use the test card if it is stored near strong electromagnetic fields.

38. Do not use the test card if it is stored near strong electric fields.

39. Do not use the test card if it is stored near strong magnetic fields.

40. Do not use the test card if it is stored near strong electromagnetic fields.

41. Do not use the test card if it is stored near strong electric fields.

42. Do not use the test card if it is stored near strong magnetic fields.

43. Do not use the test card if it is stored near strong electromagnetic fields.

44. Do not use the test card if it is stored near strong electric fields.

45. Do not use the test card if it is stored near strong magnetic fields.

46. Do not use the test card if it is stored near strong electromagnetic fields.

47. Do not use the test card if it is stored near strong electric fields.

48. Do not use the test card if it is stored near strong magnetic fields.

49. Do not use the test card if it is stored near strong electromagnetic fields.

50. Do not use the test card if it is stored near strong electric fields.

51. Do not use the test card if it is stored near strong magnetic fields.

52. Do not use the test card if it is stored near strong electromagnetic fields.

53. Do not use the test card if it is stored near strong electric fields.

54. Do not use the test card if it is stored near strong magnetic fields.

55. Do not use the test card if it is stored near strong electromagnetic fields.

56. Do not use the test card if it is stored near strong electric fields.

57. Do not use the test card if it is stored near strong magnetic fields.

58. Do not use the test card if it is stored near strong electromagnetic fields.

59. Do not use the test card if it is stored near strong electric fields.

60. Do not use the test card if it is stored near strong magnetic fields.

61. Do not use the test card if it is stored near strong electromagnetic fields.

62. Do not use the test card if it is stored near strong electric fields.

63. Do not use the test card if it is stored near strong magnetic fields.

64. Do not use the test card if it is stored near strong electromagnetic fields.

65. Do not use the test card if it is stored near strong electric fields.

66. Do not use the test card if it is stored near strong magnetic fields.

67. Do not use the test card if it is stored near strong electromagnetic fields.

68. Do not use the test card if it is stored near strong electric fields.

69. Do not use the test card if it is stored near strong magnetic fields.

70. Do not use the test card if it is stored near strong electromagnetic fields.

71. Do not use the test card if it is stored near strong electric fields.

72. Do not use the test card if it is stored near strong magnetic fields.

73. Do not use the test card if it is stored near strong electromagnetic fields.

74. Do not use the test card if it is stored near strong electric fields.

75. Do not use the test card if it is stored

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